

Rabbit Anti-Cofilin Monoclonal Antibody (CBP1131)

Cat. No.: NAB-0720-Z4274

PRODUCT OVERVIEW

Immunogen	Synthetic peptide corresponding to residues in Human Cofilin (UniProt ID: P23528).
Species Reactivity	Mouse; Rat; Human
Clonality	Monoclonal
Host Species	Rabbit
Isotype	IgG
Clone Number	CBP1131
Applications	WB; IHC-P
Research Areas	Neural Cell
Marker	Axon
Conjugation	Unconjugated

PRODUCT PROPERTIES

Formulation	Preservative: 0.01% Sodium azide Constituents: 9% PBS, 40% Glycerol, 0.05% BSA, 50% Tissue culture supernatant
Purification	Tissue culture supernatant
Purity	>95%
Shipping	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
Storage	12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution. Notes: avoid repeated freeze-thaw cycles.
Research Use Only	For research use only, not for diagnostic or therapeutic use.

TARGET

Target	Cofilin
Official Name	Cofilin
Full Name	cofilin
Alternative Names	Cofilin 1; Cofilin; Non-Muscle Isoform; Cofilin 1 (Non-Muscle); Phosphoprotein; Cofilin-1
Gene ID	1072 (Human); 12361 (Mouse); 29271 (Rat)
Uniprot ID	P23528 (Human); P18760 (Mouse); P45592 (Rat)