

iNeuMab™ Rabbit Anti-Tubulin Monoclonal Antibody (CBP1136)

Cat. No.: NAB-0720-Z4300

PRODUCT OVERVIEW

Immunogen	Synthetic peptide Within Human Tau (phospho S214). The exact sequence is proprietary. Database link: P10636
Species Reactivity	Mouse; Rat; Human
Clonality	Monoclonal
Host Species	Rabbit
Isotype	IgG
Clone Number	CBP1136
Applications	DB; WB; IHC-P
Research Areas	Neural Cell
Marker	Axon
Conjugation	Unconjugated; APC; PE; HRP; Biotin; FITC

PRODUCT PROPERTIES

Form	Liquid
Formulation	Constituent: PBS
Purification	Protein A purified
Purity	>95%
Endotoxin Level	Low Endotoxin
Shipping	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at +4°C. Avoid freeze / thaw cycle.
Storage	12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution. Notes: avoid repeated freeze-thaw cycles.
Research Use Only	NAB-0720-Z4300 is designed for research use only, not for diagnostic or therapeutic use.

TARGET

Target	Tubulin
Official Name	Tubulin
Full Name	Tubulin
Alternative Names	Tubulin