

iNeuMab™ Rat Anti-LST1 Monoclonal Antibody (CBP1198)

Cat. No.: NAB-0720-Z4564

PRODUCT OVERVIEW

Immunogen	Recombinant fragment. This information is considered to be commercially sensitive.
Species Reactivity	Mouse; Rat; Human
Clonality	Monoclonal
Host Species	Rat
Isotype	IgG2b
Clone Number	CBP1198
Applications	WB; IHC-P
Research Areas	Neural Cell
Marker	Dendrite
Conjugation	Unconjugated; APC; PE; HRP; Biotin; FITC

PRODUCT PROPERTIES

Form	Liquid
Formulation	Preservative: 0.01% Sodium azide Constituents: PBS, 0.05% BSA, 40% Glycerol
Purification	Protein A purified
Purity	>95%
Endotoxin Level	Low Endotoxin
Shipping	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage	12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution. Notes: avoid repeated freeze-thaw cycles.
Research Use Only	NAB-0720-Z4564 is designed for research use only, not for diagnostic or therapeutic use.

TARGET

Target	LST1
Official Name	LST1
Full Name	Leukocyte-specific transcript 1 protein
Alternative Names	LST1; B144; D6S49E; LST-1; leukocyte specific transcript 1
Gene ID	7940 (H uman); 16988 (Mouse); 64569 (Rat)

Uniprot ID

[O00453](#) (Human); [O08843](#) (Mouse)
