

## iNeuMab™ Rabbit Anti-SIGLEC7 Monoclonal Antibody (CBP1199)

Cat. No.: NAB-0720-Z4565

### PRODUCT OVERVIEW

<b>Immunogen</b>	Recombinant fragment Within Human Siglec 7 aa 1-400. The exact sequence is proprietary. Database link: Q9Y286
<b>Species Reactivity</b>	Human
<b>Clonality</b>	Monoclonal
<b>Host Species</b>	Rabbit
<b>Isotype</b>	IgG
<b>Clone Number</b>	CBP1199
<b>Applications</b>	FC
<b>Research Areas</b>	Neural Cell
<b>Marker</b>	Dendrite
<b>Conjugation</b>	Unconjugated; APC; PE; HRP; Biotin; FITC

### PRODUCT PROPERTIES

<b>Form</b>	Liquid
<b>Formulation</b>	Preservative: 0.01% Sodium azide Constituents: PBS, 40% Glycerol, 0.05% BSA
<b>Purification</b>	Protein A purified
<b>Purity</b>	>95%
<b>Endotoxin Level</b>	Low Endotoxin
<b>Shipping</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage</b>	12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution. Notes: avoid repeated freeze-thaw cycles.
<b>Research Use Only</b>	NAB-0720-Z4565 is designed for research use only, not for diagnostic or therapeutic use.

### TARGET

<b>Target</b>	Siglec 7
<b>Official Name</b>	SIGLEC7
<b>Full Name</b>	Sialic Acid Binding Ig Like Lectin 7

**Alternative Names**

Sialic Acid Binding Ig Like Lectin 7; QA79 Membrane Protein; SIGLEC-7; D-Siglec; AIRM-1; CDw328; AIRM1; P75; Sialic Acid Binding Ig-Like Lectin 19; Pseudogene; Adhesion Inhibitory Receptor Molecule 1; Siglec-7; Sialic Acid Binding Ig-Like Lectin; Pseudogene 2;

**Gene ID**

[27036](#) (Human)

**Uniprot ID**

[Q9Y286](#) (Human)

---