

Mouse Anti-SIGLEC7 Monoclonal Antibody (CBP1201)

Cat. No.: NAB-0720-Z4567

PRODUCT OVERVIEW

Immunogen	Recombinant fragment of Human Siglec 7 extracellular domain.
Species Reactivity	Human
Clonality	Monoclonal
Host Species	Mouse
Isotype	IgG2
Clone Number	CBP1201
Applications	WB; IHC-P
Research Areas	Neural Cell
Marker	Dendrite
Conjugation	Unconjugated

PRODUCT PROPERTIES

Formulation	Constituent: PBS
Purification	Protein G purified
Purity	>95%
Shipping	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage	12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution. Notes: avoid repeated freeze-thaw cycles.
Research Use Only	For research use only, not for diagnostic or therapeutic use.

TARGET

Target	Siglec 7
Official Name	SIGLEC7
Full Name	Sialic Acid Binding Ig Like Lectin 7
Alternative Names	Sialic Acid Binding Ig Like Lectin 7; QA79 Membrane Protein; SIGLEC-7; D-Siglec; AIRM-1; CDw328; AIRM1; P75; Sialic Acid Binding Ig-Like Lectin 19; Pseudogene; Adhesion Inhibitory Receptor Molecule 1; Siglec-7; Sialic Acid Binding Ig-Like Lectin; Pseudogene 2;
Gene ID	27036 (Human)
Uniprot ID	Q9Y286 (Human)