

## iNeuMab™ Rabbit Anti-KLK3 Monoclonal Antibody (CBP7299)

Cat. No.: NAB-0720-Z9964

### PRODUCT OVERVIEW

<b>Immunogen</b>	within Human Cathepsin K aa 300 to the C-terminus. The exact sequence is proprietary. Database link: P43235
<b>Species Reactivity</b>	Mouse; Rat; Human
<b>Clonality</b>	Monoclonal
<b>Host Species</b>	Rabbit
<b>Isotype</b>	IgG
<b>Clone Number</b>	CBP7299
<b>Applications</b>	WB
<b>Research Areas</b>	Cell Adhesion Protein
<b>Marker</b>	Proteases
<b>Conjugation</b>	Unconjugated

### PRODUCT PROPERTIES

<b>Formulation</b>	Preservative: 0.01% Sodium azide Constituents: 59% PBS, 40% Glycerol, 0.05% BSA
<b>Purification</b>	Protein A purified
<b>Purity</b>	>95%
<b>Endotoxin Level</b>	Low Endotoxin
<b>Shipping</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage</b>	12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution. Notes: avoid repeated freeze-thaw cycles.
<b>Research Use Only</b>	For research use only, not for diagnostic or therapeutic use.

### TARGET

<b>Target</b>	Prostate Specific Antigen
<b>Official Name</b>	KLK3
<b>Full Name</b>	gamma-seminoprotein or kallikrein-3
<b>Alternative Names</b>	KLK3; APS; KLK2A1; PSA; hK3; kallikrein related peptidase 3; Prostate Specific Antigen
<b>Gene ID</b>	<a href="#">354</a> (Hum an)

Uniprot ID

[P07288](#) (Human)

---