

iNeuMab™ Mouse Anti-CRMP Monoclonal Antibody (CBP4640)

Cat. No.: NAB-2103-P145

PRODUCT OVERVIEW

Immunogen	273-572 (h)
Species Reactivity	Mouse; Rat; Human
Clonality	Monoclonal
Host Species	Mouse
Isotype	mouse IgG1 κ
Clone Number	CBP4640
Applications	WB; IP; IF; ELISA

PRODUCT PROPERTIES

Purity	>95% by HPLC
Endotoxin Level	Low Endotoxin
Shipping	Gel Packs
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Research Use Only	For research use only

TARGET

Target	CRMP1
Official Name	CRMP1
Full Name	Collapsin Response Mediator Protein 1
Alternative Names	Collapsin Response Mediator Protein 1; Unc-33-Like Phosphoprotein 3; CRMP-1; DPYSL1; ULIP-3; DRP-1;
Gene ID	1400 (H uman); 12933 (Mouse); 25145 (Rat)
Uniprot ID	Q14194 (Human); P97427 (Mouse); Q62950 (Rat)
