

## Mouse Anti-PLP1 Monoclonal Antibody (CBP7353)

Cat. No.: NAB-BA-P030

### PRODUCT OVERVIEW

<b>Immunogen</b>	KLH conjugated synthetic peptide derived from rat PLP1
<b>Species Reactivity</b>	Human; Rat
<b>Clonality</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Isotype</b>	IgM
<b>Clone Number</b>	CBP7353
<b>Applications</b>	WB
<b>Research Areas</b>	Neurobiology
<b>Conjugation</b>	Unconjugated

### PRODUCT PROPERTIES

<b>Form</b>	Liquid
<b>Formulation</b>	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
<b>Concentration</b>	1mg/ml
<b>Purification</b>	affinity purified by Protein AGL
<b>Shipping</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Research Use Only</b>	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

### TARGET

<b>Target</b>	PLP1
<b>Alternative Names</b>	Proteolipid Protein 1; Lipophilin; PLP; Spastic Paraplegia 2, Uncomplicated; Major Myelin Proteolipid Protein; Pelizaeus-Merzbacher Disease; Myelin Proteolipid Protein; PLP/DM20; GPM6C; HLD1; MMPL; SPG2; PMD; MYPR_RAT;