

Mouse Anti-Ankyrin (CBP1986)

Cat. No.: NAB2007FY576

PRODUCT OVERVIEW

Immunogen	Amino acid 1-1881 of human ankyrin R fusion protein (full length)
Species Reactivity	Human
Clonality	Monoclonal
Host Species	Mouse
Isotype	IgG1
Clone Number	CBP1986
Applications	WB; IHC; ICC; IF
Conjugation	Unconjugated

PRODUCT PROPERTIES

Form	Liquid
Formulation	PBS pH7.4, 50% glycerol, 0.1% Sodium azide
Preservatives	Yes
Concentration	1 mg/ml
Purification	Protein G Purified
Shipping	The product is shipped at 4°C. Upon receipt, store it immediately at the temperature recommended below.
Storage	-20°C
Research Use Only	For research use only, not for diagnostic or therapeutic use.

TARGET

Target	ANK1
Official Name	ANK1
Alternative Names	Ankyrin-1 Antibody; Erythrocyte Ankyrin Antibody; ANK-1 Antibody; ANK Antibody; rCG_43073 Antibody
Gene ID	57463 (Human); 229715 (Mouse); 295365 (Rat)
Uniprot ID	Q86WK6 (Human); Q80ZD8 (Mouse); Q80ZD7 (Rat)