

Mouse Anti-Synaptotagmin-7 Monoclonal Antibody (CBP2085)

Cat. No.: NAB2007FY740

PRODUCT OVERVIEW

Immunogen	Amino acid 150-239 of mouse Synaptotagmin-7 fusion protein (cytoplasmic C2A domain)
Species Reactivity	Human; Mouse; Rat
Clonality	Monoclonal
Host Species	Mouse
Isotype	IgG2b
Clone Number	CBP2085
Applications	WB; IHC; ICC; IF
Conjugation	Unconjugated

PRODUCT PROPERTIES

Form	Liquid
Formulation	PBS pH7.4, 50% glycerol, 150mM NaCl, 0.02% Sodium azide
Preservatives	Yes
Concentration	1 mg/ml
Purification	Protein G Purified
Shipping	The product is shipped at 4°C. Upon receipt, store it immediately at the temperature recommended below.
Storage	-20°C
Research Use Only	For research use only, not for diagnostic or therapeutic use.

TARGET

Target	Synaptotagmin 7
Official Name	Syt7
Full Name	Synaptotagmin VII
Alternative Names	IPCA-7 Antibody, PCANAP7 Antibody, Prostate cancer-associated protein 7 Antibody, SYT-7 Antibody, SYT7 Antibody, Synaptotagmin 7 Antibody, Synaptotagmin VII Antibody, SytVII Antibody; SYT7; IPCA-7; IPCA7; PCANAP7; SYT-VII; SYTVII; synaptotagmin 7
Gene ID	9066 (H uman); 54525 (Mouse); 59267 (Rat)
Uniprot ID	O43581 (Human); Q9R0N7 (Mouse)