

Mouse Anti-DCTN2 Monoclonal Antibody (CBP2652)

Cat. No.: NAB20101726CR

PRODUCT OVERVIEW

Expression Host	Mammalian cells
Immunogen	Dynamitin (1-401 AA) fusion protein Ag1902
Species Reactivity	Human; Mouse
Clonality	Monoclonal
Host Species	Mouse
Isotype	IgG3
Clone Number	CBP2652
Applications	IHC; WB; ELISA
Research Areas	Metabolism
Conjugation	Unconjugated

PRODUCT PROPERTIES

Form	Liquid
Formulation	PBS, 0.1% sodium azide and 50% glycerol pH 7.3.
Preservatives	Yes
Concentration	1500 µg/ml
Purification	Antigen affinity purification
Purity	>95%, as determined by SDS-PAGE
Shipping	The product is shipped at 4°C. Upon receipt, store it immediately at the temperature recommended below.
Storage	Store at -20°C. Stable for one year after shipment.
Research Use Only	For research use only, not for diagnostic or therapeutic use.

TARGET

Target	DCTN2
Official Name	DCTN2
Full Name	Dynamitin 2
Alternative Names	DCTN2; dynactin 2; dynactin 2 (p50); Dynactin subunit 2; Dynamitin
Gene ID	10540 (Human); 69654 (Mouse); 299850 (Rat)
Uniprot ID	Q13561 (Human); Q99KJ8 (Mouse); Q6AYH5 (Rat)