

## Anti-CHRNb3 Monoclonal Antibody (CBP8008)

Cat. No.: NRP-0322-P1850

The antibody is a mouse monoclonal antibody raised against CHRNb3. It has been selected for its ability to recognize CHRNb3 in immunohistochemical staining and western blotting.

---

### PRODUCT OVERVIEW

<b>Species Reactivity</b>	Human
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG
<b>Clone Number</b>	CBP8008
<b>Applications</b>	WB; IHC; ICC; IP
<b>Conjugation</b>	Unconjugated

---

### PRODUCT PROPERTIES

<b>Form</b>	Liquid
<b>Formulation</b>	PBS, pH7.4, containing 0.02% NaN <sub>3</sub> , 50% glycerol.
<b>Concentration</b>	1mg/mL
<b>Purification</b>	Antigen-specific affinity chromatography followed by Protein A affinity chromatography
<b>Shipping</b>	Refrigerated packaging
<b>Storage</b>	For long term storage, the product should be stored at -20°C or lower.
<b>Research Use Only</b>	For reseaech use only, not for human or dignostic use.

---

### TARGET

<b>Target</b>	CHRNb3
<b>Official Name</b>	CHRNb3
<b>Full Name</b>	Cholinergic Receptor, Nicotinic, Beta 3
<b>Alternative Names</b>	Cholinergic Receptor Nicotinic Beta 3 Subunit; Acetylcholine Receptor; Nicotinic; Beta 3 (Neuronal); Cholinergic Receptor; Nicotinic; Beta Polypeptide 3; Cholinergic Receptor; Nicotinic; Beta 3 (Neuronal); Cholinergic Receptor; Nicotinic Beta 3;
<b>Uniprot ID</b>	Q05901