

iNeuMab™ Anti-Nicotinic Acetylcholine R alpha 4/CHRNA4 Antibody(CBP9987)

Cat. No.: NRZP-0822-ZP1933

PRODUCT OVERVIEW

Immunogen	Purified recombinant fragment of human CHRNA4 (AA: extra 29-242) expressed in E. Coli.
Species Reactivity	Human
Clonality	Monoclonal
Host Species	Mouse
Isotype	IgG1
Clone Number	CBP9987
Applications	WB; ELISA
Conjugation	Unconjugated

PRODUCT PROPERTIES

Purification	Protein A or G purified
Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	CHRNA4
Official Name	CHRNA4
Full Name	cholinergic receptor nicotinic alpha 4 subunit
Alternative Names	EBN; BFNC; EBN1; NACHR; NACRA4; NACHRA4
Gene ID	1137 (H uman)
Uniprot ID	P43681 (Human)
