

## iNeuMab™ Anti-ERR alpha/NR3B1 Antibody(CBP10092)

Cat. No.: NRZP-0822-ZP2107

### PRODUCT OVERVIEW

<b>Immunogen</b>	Purified recombinant fragment of human ESRRA (AA: 198-376) expressed in E. Coli.
<b>Species Reactivity</b>	Human; Rat
<b>Clonality</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Isotype</b>	IgG1
<b>Clone Number</b>	CBP10092
<b>Applications</b>	WB; ELISA; FC
<b>Conjugation</b>	Unconjugated

---

### PRODUCT PROPERTIES

<b>Purification</b>	Protein A or G purified
<b>Storage</b>	Store at -20°C. Do not aliquot the antibody.
<b>Research Use Only</b>	For research use only

---

### TARGET

<b>Target</b>	ESRRA
<b>Official Name</b>	ESRRA
<b>Full Name</b>	estrogen related receptor alpha
<b>Alternative Names</b>	ERR1; ERRA; ESRL1; NR3B1; ERRalpha
<b>Gene ID</b>	<a href="#">2101</a> (H uman)
<b>Uniprot ID</b>	<a href="#">P11474</a> (Human)

---