

iNeuMab™ Anti-CRMP1 Antibody (CBP11374)

Cat. No.: NRZP-0822-ZP4275

PRODUCT OVERVIEW

Immunogen	His6-tagged fusion protein corresponding to amino acids 76-592 of human CRMP1.
Species Reactivity	Human
Clonality	Monoclonal
Host Species	Mouse
Isotype	IgM Kappa
Clone Number	CBP11374
Applications	WB
Conjugation	Unconjugated

PRODUCT PROPERTIES

Purification	Protein A or G purified
Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	CRMP1
Official Name	CRMP1
Full Name	collapsin response mediator protein 1
Alternative Names	DRP1; DRP-1; CRMP-1; DPYSL1; ULIP-3
Gene ID	1400 (H uman)
Uniprot ID	Q14194 (Human)
