

iNeuMab™ Anti-NeuroD2 Antibody (CBP11398)

Cat. No.: NRZP-0822-ZP4299

PRODUCT OVERVIEW

Immunogen	NEUROD2 (NP_006151.2, 266 a.a. ~ 375 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. RTHGYCAAYETLYAAAGGGGASPDYNSSEYEGPLSPPLCLNGNFSLKQDSSPDHE KSYHYSMHYSALPGSRPTGHGLVFGSSAVRGGVHSENLISYDMHLHHDGRGPMYEE
Species Reactivity	Human
Clonality	Monoclonal
Host Species	Mouse
Isotype	IgG2a Kappa
Clone Number	CBP11398
Applications	WB; ELISA
Conjugation	Unconjugated

PRODUCT PROPERTIES

Purification	Protein A or G purified
Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	NEUROD2
Official Name	NEUROD2
Full Name	neuronal differentiation 2
Alternative Names	NDRF; DEE72; EIEE72; bHLHa1
Gene ID	4761 (H uman)
Uniprot ID	Q15784 (Human)