

iNeuMab™ Anti-CYP46A1 Antibody (CBP11458)

Cat. No.: NRZP-0822-ZP4359

PRODUCT OVERVIEW

Immunogen	CYP46A1 (NP_006659, 201 a.a. ~ 300 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. TSMLLGAQKPLSQAVKLMLEGITASRNTLAKFLPGKRKQLREVRESIRFLRQVGR DWVQRRREALKRGEV PADILTQILKAEEGAQDDEGLLDNFVTF
Species Reactivity	Human
Clonality	Monoclonal
Host Species	Mouse
Isotype	IgG1 Kappa
Clone Number	CBP11458
Applications	WB; ELISA; S-ELISA
Conjugation	Unconjugated

PRODUCT PROPERTIES

Purification	Protein A or G purified
Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	CYP46A1
Official Name	CYP46A1
Full Name	cytochrome P450 family 46 subfamily A member 1
Alternative Names	CP46; CYP46
Gene ID	10858 (Human)
Uniprot ID	Q9Y6A2(Human)