

iNeuMab™ Anti-SMPD1 Antibody (CBP11607)

Cat. No.: NRZP-0822-ZP4508

PRODUCT OVERVIEW

Immunogen	SMPD1 (AAH41164, 1 a.a. ~ 364 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MPRYGASLRQSCPRSGREQGQDGTAGAPGLLWMGLALALALALALALSDSRVLWA PAEAHPLSPQGHPARLHRIVPRLRDVFGWGNLTCPICKGLFTAINLGLKKEPNVA RVGSVAIKLCNLLKIAPPAVCRSIVHLFEDDMVEVWRRSVLSPSEACGLLLGSTC GHWDFSSWNISLPTVPKPPPKPPSPPAPGAPVSRILFLDLHWDHDYLEGTDPD CADPLCCRRGSGLPASRPGAGYWGEYSKCDLPLRTLLESLLSGLGPAGPFDMVYW TGDIPAHDVWHQTRQDQLRALTTVTALVRKFLGPVPVYPAVGNHESTPVNSFPPP FIEGNHSSRWLYEAMAKAWEPWLP AEALRTLRCI
Species Reactivity	Human
Clonality	Monoclonal
Host Species	Mouse
Isotype	IgG2a Kappa
Clone Number	CBP11607
Applications	WB; ELISA
Conjugation	Unconjugated

PRODUCT PROPERTIES

Purification	Protein A or G purified
Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	SMPD1
Official Name	SMPD1
Full Name	sphingomyelin phosphodiesterase 1
Alternative Names	aSMase; EC 3.1.4.12; NPD; ASMASE; Sphingomyelin phosphodiesterase; Acid sphingomyelinase; ASM
Gene ID	6609 (H uman)
Uniprot ID	P17405 (Human)