

## iNeuMab™ Anti-DRP1 Antibody, Clone 4E11B11

Cat. No.: NRZP-0822-ZP781

### PRODUCT OVERVIEW

<b>Immunogen</b>	Partial recombinant human DRP1 (between residues 100-250) expressed in E. coli. [Uniprot: O00429]
<b>Species Reactivity</b>	Human; Mouse; Rat; Bovine; Chinese Hamster
<b>Clonality</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Isotype</b>	IgG1
<b>Clone Number</b>	4E11B11
<b>Applications</b>	WB; ELISA; IHC; IHC-P
<b>Conjugation</b>	Unconjugated

---

### PRODUCT PROPERTIES

<b>Purification</b>	Protein A or G purified
<b>Storage</b>	Store at -20°C. Do not aliquot the antibody.
<b>Research Use Only</b>	For research use only

---

### TARGET

<b>Target</b>	DNM1L
<b>Official Name</b>	DNM1L
<b>Full Name</b>	dynamain 1 like
<b>Alternative Names</b>	DLP1; DRP1; DVLP; EMPF; OPA5; EMPF1; DYMPLE; HDYNIV
<b>Gene ID</b>	<a href="#">10059</a> (Human)
<b>Uniprot ID</b>	<a href="#">O00429(Human)</a>

---