

iNeuMab™ Anti-CRMP1 Antibody, Clone JU73-30

Cat. No.: NRZP-0822-ZP961

PRODUCT OVERVIEW

Immunogen	Recombinant protein
Species Reactivity	Human; Mouse; Rat
Clonality	Monoclonal
Host Species	Rabbit
Isotype	IgG
Clone Number	JU73-30
Applications	WB; FC; ICC; IF; IHC; IHC-P; IP
Conjugation	Unconjugated

PRODUCT PROPERTIES

Purification	Protein A or G purified
Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	CRMP1
Official Name	CRMP1
Full Name	collapsin response mediator protein 1
Alternative Names	DRP1; DRP-1; CRMP-1; DPYSL1; ULIP-3
Gene ID	1400 (H uman)
Uniprot ID	Q14194 (Human)
