

iNeuMab™ Anti-Phospho-c-Abl (CBP12753)

Cat. No.: NRZP-0922-ZP1442

The iNeuMab™ anti-Phospho-c-Abl (Tyr245) antibody detects endogenous levels of c-Abl only when phosphorylated at tyrosine 245. The antibody does not cross-react with other activated protein tyrosine kinases.

PRODUCT OVERVIEW

Expression Host	Monoclonal antibody is produced by immunizing animals with a synthetic phosphopeptide corresponding to residues surrounding Tyr245 of human c-Abl.
Species Reactivity	Human
Clonality	Monoclonal
Host Species	Rabbit
Isotype	Rabbit
Clone Number	CBP12753
Applications	WB

PRODUCT PROPERTIES

Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	ABL1
Official Name	ABL1
Full Name	Phospho-c-Abl (Tyr245)
Alternative Names	ABL; JTK7; p150; c-ABL; v-abl; CHDSKM; c-ABL1; BCR-ABL; bcr/abl
Gene ID	25 (Human)
Uniprot ID	P00519 (Human)