

iNeuMab™ Anti-Adora2a/Adenosine Receptor A2a Antibody (CBP12782)

Cat. No.: NRZP-0922-ZP1471

The iNeuMab™ anti-Adora2a/Adenosine Receptor A2a antibody (Rodent Specific) recognizes endogenous levels of total rodent A2AR protein. This antibody does not cross-react with human A2AR protein.

PRODUCT OVERVIEW

Expression Host	Monoclonal antibody is produced by immunizing animals with a synthetic peptide corresponding to residues surrounding Ala401 of mouse A2AR protein.
Species Reactivity	Mouse; Rat
Clonality	Monoclonal
Host Species	Rabbit
Isotype	Rabbit IgG
Clone Number	CBP12782
Applications	IHC; IF

PRODUCT PROPERTIES

Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	ADORA2A
Official Name	ADORA2A
Full Name	Adora2a/Adenosine Receptor A2a
Alternative Names	A2aR; RDC8; ADORA2
Gene ID	135 (Human)
Uniprot ID	P29274 (Human)