

iNeuMab™ Anti-ALDH1A1 Antibody (CBP12836)

Cat. No.: NRZP-0922-ZP1525

The iNeuMab™ anti-ALDH1A1 antibody recognizes endogenous levels of total ALDH1A1 protein.

PRODUCT OVERVIEW

Expression Host	Monoclonal antibody is produced by immunizing animals with a synthetic peptide corresponding to residues near the carboxy terminus of human ALDH1A1 protein.
Species Reactivity	Human
Clonality	Monoclonal
Host Species	Rabbit
Isotype	Rabbit IgG
Clone Number	CBP12836
Applications	WB; IP; IHC; FC

PRODUCT PROPERTIES

Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	ALDH1A1
Official Name	ALDH1A1
Full Name	ALDH1A1
Alternative Names	ALDC; ALDH1; HEL-9; HEL12; PUMB1; ALDH11; RALDH1; ALDH-E1; HEL-S-53e
Gene ID	216 (Human)
Uniprot ID	P00352 (Human)
