

iNeuMab™ Anti-CR1/CD35 Antibody (CBP13461)

Cat. No.: NRZP-0922-ZP2150

The iNeuMab™ anti-CR1/CD35 antibody recognizes endogenous levels of total CR1/CD35 protein.

PRODUCT OVERVIEW

| | |
|---------------------------|--|
| Expression Host | Monoclonal antibody is produced by immunizing animals with recombinant protein specific to the extracellular domain of human CR1/CD35 protein. |
| Species Reactivity | Human |
| Clonality | Monoclonal |
| Host Species | Rabbit |
| Isotype | Rabbit IgG |
| Clone Number | CBP13461 |
| Applications | WB |

PRODUCT PROPERTIES

| | |
|--------------------------|--|
| Storage | Store at -20°C. Do not aliquot the antibody. |
| Research Use Only | For research use only |

TARGET

| | |
|--------------------------|--------------------------------|
| Target | CR1 |
| Official Name | CR1 |
| Full Name | CR1/CD35 |
| Alternative Names | KN; C3BR; C4BR; CD35 |
| Gene ID | 1378 (H uman) |
| Uniprot ID | P17927 (Human) |
