

iNeuMab™ Anti-DRP1 Antibody (CBP13600)

Cat. No.: NRZP-0922-ZP2289

The iNeuMab™ anti-DRP1 antibody recognizes endogenous levels of total DRP1 protein.

PRODUCT OVERVIEW

Expression Host	Monoclonal antibody is produced by immunizing animals with a synthetic peptide corresponding to residues near the amino terminus of human DRP1 protein.
Species Reactivity	Human; Mouse; Rat; Monkey
Clonality	Monoclonal
Host Species	Rabbit
Isotype	Rabbit IgG
Clone Number	CBP13600
Applications	WB; IP; IF

PRODUCT PROPERTIES

Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	DNM1L
Official Name	DNM1L
Full Name	DRP1
Alternative Names	DLP1; DRP1; DVLP; EMPF; OPA5; EMPF1; DYMPLE; HDYNIV
Gene ID	10059 (Human)
Uniprot ID	Q00429(Human)
