

iNeuMab™ Anti-IDO Antibody (CBP14237)

Cat. No.: NRZP-0922-ZP2926

The iNeuMab™ anti-IDO antibody (BSA and Azide Free) recognizes endogenous levels of total IDO (IDO-1, INDO) protein. The antibody does not cross-react with IDO-2 (INDOL1). Some nonspecific staining of normal breast epithelium has been observed.

PRODUCT OVERVIEW

| | |
|---------------------------|---|
| Expression Host | Monoclonal antibody is produced by immunizing animals with recombinant human IDO protein. |
| Species Reactivity | Human |
| Clonality | Monoclonal |
| Host Species | Rabbit |
| Isotype | Rabbit IgG |
| Clone Number | CBP14237 |
| Applications | WB; IHC; IF; FC |

PRODUCT PROPERTIES

| | |
|--------------------------|--|
| Storage | Store at -20°C. Do not aliquot the antibody. |
| Research Use Only | For research use only |

TARGET

| | |
|--------------------------|--------------------------------|
| Target | IDO |
| Official Name | IDO |
| Full Name | IDO |
| Alternative Names | IDO; INDO; IDO-1 |
| Gene ID | 3620 (H uman) |
| Uniprot ID | P14902 (Human) |