

iNeuMab™ Anti-IDO Antibody (CBP14238)

Cat. No.: NRZP-0922-ZP2927

The iNeuMab™ anti-IDO antibody recognizes endogenous levels of total IDO protein. This antibody does not cross-react with mouse IDO2.

PRODUCT OVERVIEW

Expression Host	Monoclonal antibody is produced by immunizing animals with recombinant mouse IDO protein.
Species Reactivity	Mouse
Clonality	Monoclonal
Host Species	Rabbit
Isotype	Rabbit IgG
Clone Number	CBP14238
Applications	WB; IP; FC

PRODUCT PROPERTIES

Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	IDO
Official Name	IDO
Full Name	IDO
Alternative Names	IDO; INDO; IDO-1
Gene ID	3620 (H uman)
Uniprot ID	P14902 (Human)
