

Anti-VGLUT 1 Antibody (CBP17551)

Cat. No.: NRZP-1022-ZP2384

VGLUT1 is a glutamate transporter in the membrane of synaptic vesicles.

PRODUCT OVERVIEW

Immunogen	Recombinant protein corresponding to residues near the carboxy terminus of rat VGLUT1 .
Epitope	Epitop: AA 542 to 560 from rat VGLUT1.
Species Reactivity	Rat; Mouse
Clonality	Monoclonal
Host Species	Mouse
Isotype	IgG2a (κ light chain)
Clone Number	CBP17551
Applications	WB; IP; ICC; IHC; IHC-P; ELISA

PRODUCT PROPERTIES

Storage	Store at -20°C. Do not aliquot the antibody.
Research Use Only	For research use only

TARGET

Target	VGLUT1
Official Name	SLC17A7
Full Name	Vesicular glutamate transporter 1
Alternative Names	SLC17A7; BNPI; VGLUT1; solute carrier family 17 member 7
